




Realistic



Innovative




Simplicity

AFRICAN HORIZON TECHNOLOGIES Pty Ltd *developed*
Hydraspin Lab-TOC *and adapted for field use.*



The Hydraspin Lab-TOC unit is an instrument designed for the measurement of the dissolved and finely suspended organic carbon contents of water. The analysis is based on photo oxidation of organic carbon.


That is a process of ultra-violet/peroxodisulphate oxidation of the organics in the sample followed by infrared detection of the carbon dioxide liberated.



The Hydraspin Lab-TOC unit is an instrument designed for the measurement of the dissolved and finely suspended organic carbon contents of water. The analysis is based on photo oxidation of organic carbon.

That is a process of ultra-violet/peroxodisulphate oxidation of the organics in the sample followed by infrared detection of the carbon dioxide liberated.

Organic Carbon Analyser



TYPICAL APPLICATIONS



Mining



Industrial/municipal



Food & Beverage



Manufacturing



Environmental Labs



Conventional Oil & Gas

SPECIFICATIONS

Range: Linearized from 5 to 50 mg/l as carbon (dissolved organic carbon).

Response time:

Six minutes from sample introduction to result printout.

Calibration:

Pre-programmed calibration at the start of each run utilizing freshly prepared calibration standards.

Detector:

Non dispersive infra red, sensitive CO₂ gas liberated.

Sample Volume:

1 to 8 ml (Adjustable).

Outputs:

4-20 mA

Outputs:4 - 20 mA and 0- 5Volt.

Indicators:

Local LCD display

Air Supply:

Instrument quality, CO₂ Free Air, 350 ml/min at 150 Kpa.

Reagents:

Phosphoric acid solution. 05% v/v.

Potassium peroxodisulphate solution 4,0% v/v.

Reagent consumption depends on the analysis range selected but normally the acid will be consumed at a rate of 15 ml/min and the peroxodisulphate at 1 ml/min.

Temperature:

Sample 50°C maximum, instrument 40°C maximum ambient.

Power Supply:

220V 50Hz 2A.

PRINCIPLE OF OPERATION

Sample from sampling well to glass flow manifold.

Sample is mixed with Orthophosphoric acid, which reacts and removes inorganic carbon compounds.

Inorganic carbon free sample is then mixed with potassium peroxodisulphate, which aids the oxidation of organic carbon to CO₂ using ultraviolet activation .

The air/ CO₂ mixture is dried and flows through the IR detector to detect CO₂ concentration.

A digital or 4-20 mA output signal, proportional to the CO₂ concentration is fed via the Control Circuitry to the computer or server.

S.A legislation must be revised for the accurate detection of oil & grease

CONTACT US ON.....

JACQUES STEYN
+27 (0)82 781 9870
086 999 0228

jsteyn@ahtech.co.za

Excellence is not an act but a habit